



Customer name and address	
Date of contact	
Number of SUB's	Check website
Time of delivery, week from order	Check website
Version: 12	Now also micro-sparger P40
Order number	

Product name: CellMembra-300-DASbox



Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	1

Head plate layout:



Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	2

Port description:

Port #	Label	Purpose	Tubing
1	Port for DO sensor	DO.	N/A.
2	Non-invasive well	for Temp sensor.	N/A.
3	ID2,4 Dip-tube, close to bottom. - Sampling port	Sampling port.	Keep original 70 mm silicone hose (OD¼"), with swabable valve.
4	Mount condenser	DasGip hose cooler.	N/A.
5	ID6.4. - Exhaust	Exhaust line.	Keep original silicone hose, straight through connector (6060), ID6.4 (1/4") x OD11.1 (7/16") x 150 mm C-Flex hose, with OD65 (26cm ²) sterile filter.
6 O – (LA O)	ID2 Drip-guide. - Antifoam & - Perfusion -Liquid Addition overlay	Antifoam and Perfusion - Liquid Addition Overlay.	Original 70 mm silicone hose (OD¼"), Y-connection (Y230): - ID3.2 (1/8") x OD6.4 (1/4") x 500 mm C-flex hose, clamp and closed with plug. (PIP230) - ID1.6 x OD4.8 x 300 mm silicone hose, male Luer lock (MTLL220), pump tubing from EP (EP pn: 78510118), male Luer lock (MTLL230), ID3.2 (1/8") x OD6.4 (1/4") x 1000 mm C-flex hose, clamp and closed with plug. (PIP230)
6 S - (LA S)	ID2.7 Dip-tube, close to bottom. - Inoculation addition & - Perfusion - Liquid Addition	Inoculation addition & Perfusion - Liquid Addition.	Small OD10 x ID6 x 15 mm silicone hose for extra security. ID5 x OD8 x 80 mm silicone hose, female Luer-lock (FTLL240), OD25 sterile filter and male Luer plug. (MTTLP)
7	PG13.5 port	N/A.	Spare PG 13.5 port, for any 120 mm sensor, original PG13.5 plug and O-ring.
8	ID2.7 Dip-tube, set at bottom. - Harvest bottom	Bleed removal.	ID3 x OD6 x 40 mm silicone hose, straight through connector (N013-007), ID1.6 x OD4.8 x 550 silicone hose, male Luer-lock (MTLL220), pump tubing from EP (EP pn: 78510118), male Luer-lock (MTLL230), ID3.2 (1/8") x OD6.4 (1/4") x 1000 mm C-flex hose, clamp and closed with plug. (PIP230)
9 O	ID5 Dip-tube, close to bottom. - TFF Suction	TFF Suction.	ID5 x OD8 x 160 mm silicone hose to pump inlet.

Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	3

**Product code: CellMembra-300-DASbox (Part No: 30-0031-80)****PerfuseCell**

Revision 1.2

Port #	Label	Purpose	Tubing
9 S	Gassing submerged - 10 mm from the bottom.	Micro Sparger Por 3. (16-40 µm)	ID3 x OD6 x 70 mm silicone hose, Y-connection (Y230): - ID3.2 (1/8") x OD6.4 (1/4") x 500 mm C- flex hose, clamp and closed with plug. (PIP230) - ID1.6 x OD4.8 x 300 mm silicone hose, male Luer-lock (MTLL220), pump tubing from EP (EP pn: 78510118), male Luer- lock (MTLL230), ID3.2 (1/8") x OD6.4 (1/4") x 1000 mm C-flex hose, clamp and closed with plug. (PIP230)
10	Magnetic coupling	Servo motor.	N/A.
11	PG13,5 port	Selectable sensor, pH, K8.	Hamilton K8 sensor, 120 mm long.
12 A	SS level-sensor	Liquid monitoring (Estimated 65 mm from bottom, at 200ml mark)	For banana plug.
12 B	SS level-sensor	Liquid monitoring (Estimated 9 mm from bottom)	For banana plug.
13	ID5 Dip-tube, 40 mm from bottom. - TFF Return	TFF Return.	ID5 x OD8 x 280 mm silicone hose to filter barb.


Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	4

The TFF method assembly:

TFF: Custom Meissner SeptraPor®: "XFCM20C012-XX1X" - 0.2µm pore size with ID1.0 mm fibers (# lumen/fibers = 13), Length (L) = 300 mm, with a surface area of 110 cm². 1x Permeate port.

Permeate; OD8 x ID5 x 40 mm silicone hose, Y-connection (260-6005):



- A. **Wash-Line:** OD8 x ID5 x 100 mm silicone hose, clamp, female Luer-Lock (FTLL360) and swabable valve.
- B. **Harvest-Line:** OD8 x ID5 x 100 mm silicone hose, clamp, straight through reducer (3050), ID1.6 x OD4.8 x 450 mm silicone hose, male Luer-Lock (MTLL220), Bioprene pump-head tubing (ID1.0 mm x OD3.1 mm), male Luer-Lock (MTLL230), ID3.2 (1/8") x OD6.4 (1/4") x 600 mm C-flex hose, clamp, Y-connection (240-6005):
 - 1. ID3.2 (1/8") x OD6.4 (1/4") x 100 mm C-flex hose, clamp, female Luer-Lock (FTLL230), and swabable valve.
 - 2. ID3.2 (1/8") x OD6.4 (1/4") x 400 mm C-flex hose, clamp and plug. (PIP230)





Mounting Port-Addition-Mount (PAM-TFF)

Insert SS-rods in the PAM-TFF 3D-print and glue in place with a small amount of UV-glue as shown in picture below. (shortest SS-rod in the left side – Port 12A)

Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	5




Position the PAM-TFF so that both dip-tubes and SS-rods are inside the vessel. Then dispense a suitable amount of silicone on to the headplate filling up the area of the PAM-TFF. Carefully press down the PAM-TFF until its “legs” touches the vessel sealing it in place.




Remove any excess silicone if needed and attach the Clio vessel support. Make sure the snapfits lock in place underneath the vessel lid and let dry for 24 hours.

Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	6




Clio SUP PC tube support:




Attached OD10 x 240 mm PC tube for support.

Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	7

Hollow-Fiber-Filter orientation and green "PerfuseCell Cap":



The HFF to be mounted with the permeate barb parallel to the assembly pointing slightly to the right of the support PC tube, as shown in the picture above:



When assembly is complete slide on the green "PerfuseCell Cap".

Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	8


Single-Use-Pump:

The CM60 SUP used from 2017 up to 2025 was a machined Polycarbonate body with epoxy glued on dome and fluid velocity limited by the Spectrum HFF Luer-Lock connections.

From 2025 (ver 10 onwards) the Nylon medical grade Clio M80-30 SUP integrate the SeptraPor® “XFCM20C012-XX1X” Hollow-Fiber-Filter cartridge directly in the pump body eliminating Luer-Lok for reduced shear stress. Medical grade molded silicone diaphragm.

SeptraPor® HFF spec is: 1 mm lumen, 13 lumen, 0.2 µm / 200 nm pore, 110 cm² area.

<https://perfusecell.com/perfusion-bioreactors/cellmembra-tff-integrated-30/cellmembra-300-bioblu>



Level mark:

Level mark needs to be placed in front of “Exhaust”.

Strips:

Strips must be applied to all hose barbs.

All strips attached to the headplate to be positioned so that the square lock faces outward, away from the center of the headplate and toward the outer diameter of the vessel.

Designed and manufactured by

PerfuseCell-Cronus-PCS A/S

Noerrelundvej 10

DK-2730 Herlev

Denmark

www.perfusecell.com

Dokumentnummer:	Emne:				
K08.2.3.1.7	Production ordre formular – standard product				
Version:	Redaktør:	Dokumentansvarlig:	Revised by:	Godkendt (dato):	Side:
1	PS	PS	PS	2025-03-31	9